

FIT3077 Team Indomie A2

System Manual

Follow these instructions step by step to use our system.

1. Clone the system repository in your machine.
2. Open the file `PracIdView.java`
3. Run the contents of the file and a window should popup.
4. Enter your practitioner identifier and click on enter.
5. Please wait for the system to complete the fetch of your patients. A new window of your patients will appear once completed.
6. To add a patient to monitor, select their row on the patient list. Click "Add Patient to Monitor". Please note that a patient can only be monitored if they have a cholesterol reading.
7. Click on "Cholesterol Level" to open the cholesterol monitor.
8. To view the patients personal information, click on the row and their data will appear on the text fields to the right.
9. To add more patients to monitor, follow step 6.
10. The cholesterol level of your patients are highlighted red if it is higher than the table's cholesterol average.
11. To set the monitor's update interval time, enter an integer in the "Set Timer" text field and click on "Set Timer" to update.
12. To remove a patient from the monitor, click on the row of the patient to remove. Check on the patient's data field that it is the correct patient, and click "Remove Patient" to proceed.
13. To close the monitor, simply close the window of the cholesterol monitor.
14. To stop the system, close the window of your patient list.